


Product datasheet

Anti-ZIP-12 antibody ab64463

画像数 1

製品の概要

製品名	Anti-ZIP-12 antibody
製品の詳細	Rabbit polyclonal to ZIP-12
由来種	Rabbit
アプリケーション	適用あり: WB, ELISA
種交差性	交差種: Human 交差が予測される動物種: Horse 
免疫原	Synthetic peptide (Human)(Peptide available as ab113082 .)
ポジティブ・コントロール	Jurkat cells

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C.
バッファー	Preservative: 0.09% Sodium azide Constituents: 2% Sucrose, PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab64463** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 0.5 µg/ml. Detects a band of approximately 68 kDa (predicted molecular weight: 73 kDa). Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.


アプリ ケーション	Abreviews	特記事項
--------------	-----------	------

ELISA
Use at an assay dependent concentration.
Titre using peptide based assay 1:62500.

ターゲット情報

機能	Acts as a zinc-influx transporter (Potential). May be partly involved in the outbreak of schizophrenia.
組織特異性	Expressed in brain and eye.
配列類似性	Belongs to the ZIP transporter (TC 2.A.5) family.
細胞内局在	Membrane.

画像



Anti-ZIP-12 antibody (ab64463) at 0.5 µg/ml + Jurkat cell lysate at 10 µg

Secondary

HRP conjugated anti-Rabbit IgG at 1/50000 dilution

Predicted band size: 73 kDa

Observed band size: 68 kDa

[why is the actual band size different from the predicted?](#)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors