


## Product datasheet

# Anti-ZNF541 antibody ab66264

### 画像数 1

#### 製品の概要

製品名	Anti-ZNF541 antibody
製品の詳細	Rabbit polyclonal to ZNF541
由来種	Rabbit
アプリケーション	適用あり: WB, ELISA
種交差性	交差種: Human 交差が予測される動物種: Mouse, Rat, Horse, Guinea pig, Cow, Cat, Dog 
免疫原	A 14-22 amino acid region within synthetic peptide: ASPSQVAMAS FSSAGPPADP SKSKLTIFSR IQGGNYRLP HPVKEENVAG, corresponding to internal sequence amino acids 360-409 of Human ZNF541 <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
ポジティブ・コントロール	Jurkat cell lysate.

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: 0.09% Sodium azide Constituents: 2% Sucrose, PBS
精製度	Protein A purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab66264** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリ ケーション	Abreviews	特記事項
WB		Use a concentration of 1.25 µg/ml. Detects a band of approximately 50 kDa (predicted molecular weight: 87 kDa). Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.
ELISA		Use at an assay dependent concentration. Titre using peptide based assay 1/1562500.

## ターゲット情報

機能	Component of some chromatin remodeling multiprotein complex that plays a role during spermatogenesis.
配列類似性	Contains 5 C2H2-type zinc fingers. Contains 1 ELM2 domain. Contains 1 SANT domain.
細胞内局在	Nucleus.

## 画像



**All lanes :** Anti-ZNF541 antibody (ab66264) at 1.25 µg/ml

**Lane 1 :** Molecular weight markers

**Lane 2 :** Jurkat cell lysate at 10 µg

### Secondary

**All lanes :** HRP conjugated anti-Rabbit IgG at 1/50000 dilution

**Predicted band size:** 87 kDa

**Observed band size:** 50 kDa

[why is the actual band size different from the predicted?](#)

Antibody diluted in 5% skim milk / PBS buffer. Gel concentration 12%.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors