

## Product datasheet

# Anti-C8orf72 antibody ab67209

### 画像数 1

#### 製品の概要

製品名	Anti-C8orf72 antibody
製品の詳細	Mouse polyclonal to C8orf72
由来種	Mouse
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Full length protein (Human)

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: Whole serum
精製度	Whole antiserum
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab67209** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/1000. Detects a band of approximately 51 kDa (predicted molecular weight: 41 kDa).

#### ターゲット情報

機能	May be involved in tumor progression.
----	---------------------------------------

### 組織特異性

Detected in thyroid, spleen and testis, and at lower levels in stomach, spinal cord, lymph node, trachea, adrenal gland, prostate, ovary and intestine.


### 配列類似性

Belongs to the FAM110 family.

### 細胞内局在

Cytoplasm. Cytoplasm > cytoskeleton > centrosome.

### 画像



**All lanes** : Anti-C8orf72 antibody (ab67209)  
at 1/500 dilution

**Lane 1** : C8orf72 transfected 293T cell lysate

**Lane 2** : Non transfected 293T cell lysate

**Predicted band size:** 41 kDa

**Observed band size:** 51 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors