

Product datasheet

Anti-IFNA4 antibody ab67385

画像数 2

製品の概要

製品名	Anti-IFNA4 antibody
製品の詳細	Mouse polyclonal to IFNA4
由来種	Mouse
アプリケーション	適用あり: WB, ICC/IF
種交差性	交差種: Human
免疫原	Full length protein, corresponding to amino acids 1-189 of Human IFNA4
ポジティブ・コントロール	Transfected 293Tcell lysate (WB) HeLa cells (ICC/IF)

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: None Constituents: Whole serum
精製度	Whole antiserum
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab67385** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		
ICC/IF		

追加情報

ICC/IF: Use at a concentration of 10 µg/ml. Purification with Protein A is recommended prior to proceeding with the ICC/IF protocol.

WB: 1/500 - 1/1000. Detects a band of approximately 20 kDa (predicted molecular weight: 22 kDa).

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

ターゲット情報


関連性

Produced by macrophages, IFN alpha have antiviral activities. Interferon stimulates the production of two enzymes: a protein kinase and an oligoadenylate synthetase.

細胞内局在


Secreted

画像



ab67385, at a concentration 10 μg/ml, staining IFNA4 in HeLa cell by Immunofluorescence.

Immunocytochemistry/ Immunofluorescence - Anti-IFNA4 antibody (ab67385)



Western blot - Anti-IFNA4 antibody (ab67385)

All lanes : Anti-IFNA4 antibody (ab67385) at 1/500 dilution

Lane 1 : IFNA4 transfected 293T cell lysate

Lane 2 : non transfected 293T cell lysate

Lysates/proteins at 25 μg per lane.

Secondary

All lanes : Goat Anti-Mouse IgG (H&L)-HRP

Conjugate secondary antibody at 1/2500 dilution

Predicted band size: 22 kDa

Observed band size: 18,20 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors