

Product datasheet

Anti-NPM2 antibody ab67403

★★★★☆ 1 Abreviews 1 References 画像数 1

製品の概要

製品名	Anti-NPM2 antibody
製品の詳細	Mouse polyclonal to NPM2
由来種	Mouse
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Full length protein, corresponding to amino acids 1-214 of Human NPM2
ポジティブ・コントロール	Transfected 293T cells

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: None Constituents: PBS, pH 7.2
精製度	Protein A purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab67403** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/1000. Detects a band of approximately 34 kDa (predicted molecular weight: 24 kDa).

ターゲット情報

機能 Core histones chaperone involved in chromatin reprogramming, specially during fertilization and

early embryonic development. Probably involved in sperm DNA decondensation during fertilization.

配列類似性

Belongs to the nucleoplasmin family.


ドメイン

The acidic tract A2 mediates histone binding.

細胞内局在

Nucleus. Found in the oocyte nucleus before nuclear membrane breakdown, after which it is redistributed to the cytoplasm.

画像



All lanes : Anti-NPM2 antibody (ab67403) at 1/500 dilution

All lanes : NPM2 transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

Secondary

All lanes : Goat Anti-Mouse IgG (H&L)-HRP Conjugate at 1/2500 dilution

Predicted band size: 24 kDa

Observed band size: 34 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors