

Product datasheet

Anti-BYSL antibody ab67414

画像数 1

製品の概要

製品名	Anti-BYSL antibody
製品の詳細	Mouse polyclonal to BYSL
由来種	Mouse
特異性	ab67414 is reactive against mammalian transfected lysate.
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Full length protein, corresponding to amino acids 1-325 of Human BYSL
ポジティブ・コントロール	BYSL Transfected 293T cell lysate

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: 1X PBS, pH 7.5
精製度	Protein A purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション


Our [Abpromise guarantee](#) covers the use of **ab67414** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/1000. Detects a band of approximately 37 kDa (predicted molecular weight: 50 kDa).

機能	Required for processing of 20S pre-rRNA precursor and biogenesis of 40S ribosomal subunits. May be required for trophinin-dependent regulation of cell adhesion during implantation of human embryos.
組織特異性	Found in the placenta from the sixth week of pregnancy. Was localized in the cytoplasm of the syncytiotrophoblast in the chorionic villi and in endometrial decidual cells at the uteroplacental interface. After week 10, the level decreased and then disappeared from placental villi.
配列類似性	Belongs to the bystin family.
細胞内局在	Cytoplasm. Nucleus > nucleolus. Associated with 40S ribosomal subunits.

画像



Western blot - Anti-BYSL antibody (ab67414)

All lanes : Anti-BYSL antibody (ab67414) at 1/500 dilution

Lane 1 : BYSL transfected 293T lysate

Lane 2 : Non transfected 293T lysate

Lysates/proteins at 25 µg per lane.

Secondary

All lanes : Goat Anti-Mouse IgG (H&L)-HRP Conjugate at 1/2500 dilution

Predicted band size: 50 kDa

Observed band size: 37 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors