

Product datasheet

Anti-PCDHA8 antibody ab67611

画像数 1

製品の概要

製品名	Anti-PCDHA8 antibody
製品の詳細	Mouse polyclonal to PCDHA8
由来種	Mouse
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Full length protein (Human, aa 1-814)

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Constituent: Whole serum
精製度	Whole antiserum
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab67611** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/1000. Detects a band of approximately 95 kDa (predicted molecular weight: 103 , 89 kDa).

ターゲット情報


関連性 PCDHA8 is most likely to play a critical role in the establishment and function of specific cell-cell connections in the brain. Alternative splicing has been observed and additional variants have

been suggested but their full-length nature has yet to be determined. Two different isoforms are known.

細胞内局在

Cell membrane; Single-pass type I membrane protein

画像



Western blot - Anti-PCDHA8 antibody (ab67611)

All lanes : Anti-PCDHA8 antibody (ab67611) at 1/500 dilution

Lane 1 : PCDHA8 transfected 293T cell lysate

Lane 2 : Non-transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

Secondary

All lanes : Goat Anti-Mouse IgG (H&L)-HRP conjugated at 1/2500 dilution

Predicted band size: 103 , 89 kDa

Observed band size: 95 kDa

[why is the actual band size different from the predicted?](#)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors