

Product datasheet

Anti-TEKT3 antibody ab67779

画像数 1

製品の概要

製品名	Anti-TEKT3 antibody
製品の詳細	Mouse polyclonal to TEKT3
由来種	Mouse
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Full length protein (Human)

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: Whole serum
精製度	Whole antiserum
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab67779** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/1000. Detects a band of approximately 57 kDa (predicted molecular weight: 57 kDa).

ターゲット情報

機能	Structural component of ciliary and flagellar microtubules. Forms filamentous polymers in the
----	---

walls of ciliary and flagellar microtubules. Required for progressive sperm mobility.


#### 配列類似性

Belongs to the tektin family.

#### 細胞内局在

Cytoplasm > cytoskeleton > cilium axoneme. Cytoplasm > cytoskeleton > flagellum axoneme.

#### 画像



**All lanes** : Anti-TEKT3 antibody (ab67779) at 1/500 dilution

**Lane 1** : TEKT3 transfected 293T cell lysate

**Lane 2** : Non-transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

#### Secondary

**All lanes** : Goat Anti-Mouse IgG (H&L)-HRP conjugated at 1/2500 dilution

**Predicted band size:** 57 kDa

**Observed band size:** 57 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

#### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

#### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors