

Product datasheet

Anti-UBE2Z antibody ab68001

画像数 1

製品の概要

製品名	Anti-UBE2Z antibody
製品の詳細	Mouse polyclonal to UBE2Z
由来種	Mouse
特異性	This antibody is reactive against mammalian transfected lysate.
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	UBE2Z (NP_075567.1, 1 a.a. ~ 246 a.a) full-length human protein.
ポジティブ・コントロール	UBE2Z transfected 293T cell lysate

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	pH: 7.20 Constituent: PBS
精製度	Protein A purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション


Our [Abpromise guarantee](#) covers the use of **ab68001** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/1000. Detects a band of approximately 31 kDa (predicted molecular weight: 38 kDa).

機能	Catalyzes the covalent attachment of ubiquitin to other proteins (By similarity). Specific substrate for UBE1L2, not charged with ubiquitin by UBE1. May be involved in apoptosis regulation.
組織特異性	Widely expressed. Highly in placenta, pancreas, spleen and testis.
パスウェイ	Protein modification; protein ubiquitination.
配列類似性	Belongs to the ubiquitin-conjugating enzyme family.
細胞内局在	Cytoplasm. Nucleus.

画像



Western blot - Anti-UBE2Z antibody (ab68001)

All lanes : Anti-UBE2Z antibody (ab68001) at 1/500 dilution

Lane 1 : UBE2Z transfected 293T cell lysate

Lane 2 : Non-transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

Secondary

All lanes : Goat Anti-Mouse IgG (H&L)-HRP

Conjugate

at 1/2500 dilution

Predicted band size: 38 kDa

Observed band size: 31 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors