

## Product datasheet

# Anti-OTUB2 antibody ab68010

★★★★☆ 1 Abreviews [画像数 1](#)

### 製品の概要

製品名	Anti-OTUB2 antibody
製品の詳細	Mouse polyclonal to OTUB2
由来種	Mouse
特異性	This antibody is reactive against mammalian transfected lysate.
アプリケーション	<b>適用あり:</b> WB
種交差性	<b>交差種:</b> Human
免疫原	OTUB2 (NP_075601.1, 1 a.a. ~ 234 a.a) full-length human protein.
ポジティブ・コントロール	OTUB2 transfected 293T cell lysate

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: Whole serum
精製度	Whole antiserum
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション


Our [Abpromise guarantee](#) covers the use of **ab68010** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB	★★★★☆	1/500 - 1/1000. Detects a band of approximately 25 kDa (predicted molecular weight: 27 kDa).

<b>機能</b>	Hydrolase that can remove conjugated ubiquitin from proteins in vitro and may therefore play an important regulatory role at the level of protein turnover by preventing degradation. Mediates deubiquitination of both 'Lys-48'- and 'Lys-63'-linked polyubiquitin chains, with a preference for 'Lys-63'-linked polyubiquitin chains.
<b>組織特異性</b>	Widely expressed. Expressed at higher level in brain.
<b>配列類似性</b>	Belongs to the peptidase C65 family. Contains 1 OTU domain.

## 画像



**All lanes :** Anti-OTUB2 antibody (ab68010) at 1/500 dilution

**Lane 1 :** OTUB2 transfected 293T cell lysate

**Lane 2 :** Non-transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

### Secondary

**All lanes :** Goat Anti-Mouse IgG (H&L)-HRP  
Conjugate at 1/2500 dilution

**Predicted band size:** 27 kDa

**Observed band size:** 25 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors