

Product datasheet

Anti-RBM6 antibody ab68297

画像数 1

製品の概要

製品名	Anti-RBM6 antibody
製品の詳細	Mouse polyclonal to RBM6
由来種	Mouse
特異性	This antibody is reactive against mammalian transfected lysate.
アプリケーション	<b>適用あり:</b> WB
種交差性	<b>交差種:</b> Human
免疫原	RBM6 (AAH46643.1, 1 a.a. ~ 601 a.a) full-length human protein.
ポジティブ・コントロール	RBM6 transfected 293T cell lysate

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: 1X PBS, pH 7.2
精製度	Protein A purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション


Our [Abpromise guarantee](#) covers the use of **ab68297** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/1000. Detects a band of approximately 70 kDa.

<b>機能</b>	Specifically binds poly(G) RNA homopolymers in vitro.
<b>組織特異性</b>	Ubiquitous in adults.
<b>配列類似性</b>	Contains 1 G-patch domain. Contains 1 RRM (RNA recognition motif) domain.
<b>細胞内局在</b>	Nucleus.

## 画像



**All lanes :** Anti-RBM6 antibody (ab68297) at 1/500 dilution

**Lane 1 :** RBM6 transfected 293T cell lysate

**Lane 2 :** Non-transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

### Secondary

**All lanes :** Goat Anti-Mouse IgG (H&L)-HRP Conjugate at 1/2500 dilution

**Observed band size:** 70 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors