

Product datasheet

Anti-GNL2 antibody ab68783

画像数 1

製品の概要

製品名	Anti-GNL2 antibody
製品の詳細	Mouse polyclonal to GNL2
由来種	Mouse
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Full length GNL2 protein (Human) (NP_037417)
ポジティブ・コントロール	GNL2 transfected 293T cell lysate

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: Whole serum
精製度	Whole antiserum
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab68783** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/1000. Detects a band of approximately 100 kDa (predicted molecular weight: 84 kDa).

ターゲット情報

機能	GTPase that associates with pre-60S ribosomal subunits in the nucleolus and is required for their nuclear export and maturation.
組織特異性	Ubiquitously expressed at relatively low levels in all human tissues tested, with the highest level of expression in the testes.
配列類似性	Belongs to the MMR1/HSR1 GTP-binding protein family. NOG2 subfamily. Contains 1 G (guanine nucleotide-binding) domain.
細胞内局在	Nucleus > nucleolus.

画像



All lanes : Anti-GNL2 antibody (ab68783) at 1/500 dilution

Lane 1 : GNL2 transfected 293T cell lysate at 25 µg

Lane 2 : Non transfected 293T cell lysate

Secondary

All lanes : Goat Anti-Mouse IgG (H&L)-HRP Conjugate secondary antibody at 1/2500 dilution

Developed using the ECL technique.

Predicted band size: 84 kDa

Observed band size: 100 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors