

Product datasheet

Anti-Neuronal calcium-binding protein antibody ab69204

画像数 1

製品の概要

製品名	Anti-Neuronal calcium-binding protein antibody
製品の詳細	Mouse polyclonal to Neuronal calcium-binding protein
由来種	Mouse
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Full length human NECAB3 protein (AAH47673.1).
ポジティブ・コントロール	NECAB3 transfected 293T lysate.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: Whole serum
精製度	Whole antiserum
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab69204** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/1000. Detects a band of approximately 36 kDa (predicted molecular weight: 41 kDa).

ターゲット情報

機能	Inhibits the interaction of APBA2 with beta-amyloid precursor protein (APP), and hence allows formation of beta-amyloid.
組織特異性	Strongly expressed in heart and skeletal muscle, moderately in brain and pancreas.
配列類似性	Contains 1 ABM domain. Contains 1 EF-hand domain.
翻訳後修飾	Phosphorylated by NEK2.
細胞内局在	Golgi apparatus.

画像



Western blot - Anti-Neuronal calcium-binding protein antibody (ab69204)

All lanes : Anti-Neuronal calcium-binding protein antibody (ab69204) at 1/500 dilution

Lane 1 : NECAB3 transfected 293T lysate

Lane 2 : Non-transfected 293T lysate

Lysates/proteins at 25 µg per lane.

Secondary

All lanes : Goat Anti-Mouse IgG (H&L)-HRP at 1/2500 dilution

Predicted band size: 41 kDa

Observed band size: 36 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors