

## Product datasheet

# Anti-STK11IP antibody ab70441

### 画像数 1

#### 製品の概要

製品名	Anti-STK11IP antibody
製品の詳細	Mouse polyclonal to STK11IP
由来種	Mouse
特異性	This antibody is reactive against mammalian transfected lysate.
アプリケーション	<b>適用あり:</b> WB
種交差性	<b>交差種:</b> Human
免疫原	Full length protein, corresponding to amino acids 1-1099 of Human STK11IP
ポジティブ・コントロール	Transfected 293T cell lysate.

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: 1X PBS, pH 7.2
精製度	Protein G purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab70441** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
----------	-----------	------

WB

#### 追加情報

WB: 1/500 - 1/1000. Detects a band of approximately 150 kDa (predicted molecular weight: 121 kDa).

Detection Antibody.


Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

## ターゲット情報

機能	May regulate STK11 function by controlling its subcellular localization.
配列類似性	Belongs to the STK11IP family. Contains 8 LRR (leucine-rich) repeats.
細胞内局在	Cytoplasm. Some cells show granular or punctuate expression. Colocalizes with STK11 and SMAD4 in granular or punctuate structures.

## 画像



**All lanes** : Anti-STK11IP antibody (ab70441)  
at 1/500 dilution

**Lane 1** : STK11IP transfected 293T cell lysate

**Lane 2** : Non-transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

### Secondary

**All lanes** : Goat Anti-Mouse IgG (H&L)-HRP  
Conjugate at 1/2500 dilution

**Predicted band size:** 121 kDa

**Observed band size:** 150 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors