

Product datasheet

Goat Anti-Syrian Hamster IgG H&L preadsorbed ab7141

製品の概要

製品名	Goat Anti-Syrian Hamster IgG H&L preadsorbed
由来種	Goat
ターゲット生物種	Syrian hamster
アプリケーション	適用あり: IHC-P, IHC-Fr, ICC/IF, IP, Immunodiffusion, Conjugation, WB, ELISA
吸着処理血清	Chicken, Cow, Goat, Guinea pig, Horse, Human, Mouse, Rabbit, Rat, Sheep more details
免疫原	Hamster IgG whole molecule

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 7.60 Preservative: 0.01% Sodium azide Constituents: 0.88% Sodium chloride, 0.42% Potassium phosphate
精製度	Affinity purified
特記事項 (精製)	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Hamster IgG coupled to agarose beads.
ポリ/モノ	ポリクローナル

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab7141** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
IHC-P		Use at an assay dependent concentration.
IHC-Fr		Use at an assay dependent concentration.
ICC/IF		1/1000 - 1/5000.

アプリケーション	Abreviews	特記事項
IP		Use at an assay dependent concentration.
Immunodiffusion		Use at an assay dependent concentration.
Conjugation		Use at an assay dependent concentration.
WB		1/5000 - 1/50000.
ELISA		1/800000.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors