


Product datasheet

Anti-UNC93B antibody ab72123

画像数 1

製品の概要

製品名	Anti-UNC93B antibody
製品の詳細	Rabbit polyclonal to UNC93B
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Mouse, Rat, Human 交差が予測される動物種: Dog 
免疫原	A Synthetic peptide (Mouse) corresponding to residues 500-550 of mouse UNC93B. (NP_062322).
ポジティブ・コントロール	Human spleen lysate Mouse spleen lysate

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: 0.05% Sodium Azide Constituents: 0.05% BSA, PBS
精製度	Protein G purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab72123** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
WB		Use a concentration of 3 - 5 µg/ml. Detects a band of approximately 67 kDa (predicted molecular weight: 67 kDa).

ターゲット情報

関連性 UNC93B is a protein with similarity to the *C. elegans* unc93 protein. The Unc93 protein is involved in the regulation or coordination of muscle contraction in the worm. Defects in UNC93B are associated with susceptibility to herpes simplex encephalitis (HSE) which is a rare complication of human herpesvirus 1 (HHV-1) infection.

細胞内局在 Membrane; Multi pass membrane protein

画像



All lanes : Anti-UNC93B antibody (ab72123)
at 4 µg/ml

Lane 1 : Human spleen lysate

Lane 2 : Human spleen lysate plus
immunising peptide

Lane 3 : Mouse spleen lysate

Predicted band size: 67 kDa

Observed band size: 67 kDa

Additional bands at: 30 kDa. We are unsure
as to the identity of these extra bands.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors