


Product datasheet

Anti-SNX2 antibody ab72197

画像数 1

製品の概要

製品名	Anti-SNX2 antibody
製品の詳細	Rabbit polyclonal to SNX2
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Mouse, Rat, Human, Monkey 交差が予測される動物種: Chicken, Cow, Xenopus laevis 
免疫原	Synthetic peptide corresponding to a portion of amino acids 120-170 of Human SNX2
ポジティブ・コントロール	Mouse placenta lysate.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: 0.05% Sodium azide Constituents: 0.05% BSA, PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab72197** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 0.25 - 0.5 µg/ml. Detects a band of approximately 70 kDa (predicted molecular weight: 58 kDa).

ターゲット情報

機能	May be involved in several stages of intracellular trafficking. Component of the retromer complex, a complex required to retrieve lysosomal enzyme receptors (IGF2R and M6PR) from endosomes to the trans-Golgi network. Interacts with membranes containing phosphatidylinositol 3-phosphate (PtdIns(3P)) or phosphatidylinositol 3,5-bisphosphate (PtdIns(3,5)P2).
配列類似性	Belongs to the sorting nexin family. Contains 1 PX (phox homology) domain.
細胞内局在	Early endosome membrane.

画像



All lanes : Anti-SNX2 antibody (ab72197) at 0.25 µg/ml

Lane 1 : Mouse placenta lysate

Lane 2 : Mouse placenta lysate with immunising peptide

Predicted band size: 58 kDa

Observed band size: 70 kDa

[why is the actual band size different from the predicted?](#)

Additional bands at: 15 kDa, 30 kDa, 35 kDa. We are unsure as to the identity of these extra bands.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors