

Product datasheet

Anti-ERAS antibody ab72401

1 References [画像数 1](#)

製品の概要

製品名	Anti-ERAS antibody
製品の詳細	Rabbit polyclonal to ERAS
由来種	Rabbit
アプリケーション	適用あり: IHC-P
種交差性	交差種: Rat, Human 交差が予測される動物種: Mouse
免疫原	Synthetic peptide: VETSAKTRQGVEEAFSL , corresponding to amino acids 178-195 of Human ERAS <div style="text-align: right;"> Run BLAST with Run BLAST with </div>
ポジティブ・コントロール	Rat brain tissue

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at -20°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.025% Sodium Azide Constituents: 50% Glycerol, Whole serum
精製度	Whole antiserum
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab72401** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
IHC-P		1/100.

ターゲット情報

機能	Ras proteins bind GDP/GTP and possess intrinsic GTPase activity. Plays an important role in the tumor-like growth properties of embryonic stem cells.
配列類似性	Belongs to the small GTPase superfamily. Ras family.
細胞内局在	Cell membrane.

画像



Immunohistochemical analysis of ERAS in paraffin-embedded formalin-fixed rat brain tissue using an isotype control (top, left) and ab72401 (bottom left, and right) at 1:100.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-ERAS antibody (ab72401)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors