

Product datasheet

Anti-JRK antibody ab72683

画像数 1

製品の概要

製品名	Anti-JRK antibody
製品の詳細	Mouse polyclonal to JRK
由来種	Mouse
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Full length human JRK protein (AAH43351.1)
ポジティブ・コントロール	JRK transfected 293T lysate

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: PBS, pH 7.2
精製度	Protein A purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab72683** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
WB		1/500 - 1/1000. Predicted molecular weight: 57 kDa.

ターゲット情報

関連性 The JRK gene is the human homolog of the mouse jerky gene. The encoded protein has similarity to several nuclear regulatory proteins, including centromere protein B, which indicates that it might function as a DNA binding protein. Insertional inactivation of this gene in transgenic mice results in epileptic seizures. JRK is a strong candidate for Childhood Absence Epilepsy (CAE) as it has been mapped to the same chromosomal location (8q24.3). Two transcript variants encoding different isoforms have been found for this gene.

細胞内局在 Nuclear

画像



All lanes : Anti-JRK antibody (ab72683) at 1/500 dilution

Lane 1 : JRK transfected 293T lysate

Lane 2 : Non transfected 293T lysate

Lysates/proteins at 25 µg per lane.

Secondary

All lanes : Goat Anti-Mouse IgG (H&L)-HRP Conjugate at 1/2500 dilution

Predicted band size: 57 kDa

Observed band size: 62 kDa

[why is the actual band size different from the predicted?](#)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors