

## Product datasheet

# Anti-Cytokeratin 6 antibody [SPM269], prediluted ab75690

### 画像数 1

#### 製品の概要

製品名	Anti-Cytokeratin 6 antibody [SPM269], prediluted
製品の詳細	Mouse monoclonal [SPM269] to Cytokeratin 6, prediluted
由来種	Mouse
アプリケーション	適用あり: IHC-P
種交差性	交差種: Mouse, Human
免疫原	A synthetic peptide from the C terminus of human Cytokeratin 6
ポジティブ・コントロール	Human skin tissue (suprabasal cells of the outer hair follicle root sheath in skin).

#### 製品の特性

製品の状態	Prediluted
保存方法	Shipped at 4°C. Store at +4°C.
特記事項(精製)	ab75690 is purified from ascites fluid by Protein A.
ポリ/モノ	モノクローナル
クローン名	SPM269
アイソタイプ	IgG2a
軽鎖の種類	kappa

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab75690** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
IHC-P		

#### 追加情報

IHC-P: Ready to use for 30 minutes at room temperature.  
 Staining of formalin fixed tissues requires boiling tissue sections in 10mM citrate buffer, pH 6.0 for 10 min followed by cooling at room temperature for 20 minutes.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

## ターゲット情報

---

### 組織特異性

Constitutively expressed in distinct types of epithelia such as those in oral mucosa, esophagus, papillae of tongue and hair follicle outer root sheath.


### 配列類似性

Belongs to the intermediate filament family.

---

## 画像

---



ab75690, prediluted, staining Cytokeratin 6 in formalin fixed, paraffin embedded human skin tissue by Immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Cytokeratin 6 antibody [SPM269], prediluted (ab75690)

---

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors