

## Product datasheet

# Anti-FSH antibody [SPM107] ab75697

### 画像数 1

#### 製品の概要

製品名	Anti-FSH antibody [SPM107]
製品の詳細	Mouse monoclonal [SPM107] to FSH
由来種	Mouse
アプリケーション	適用あり: IHC-P
種交差性	交差種: Human
免疫原	Full length purified Human FSH
ポジティブ・コントロール	Anterior pituitary

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	pH: 7.40 Preservative: 0.1% Sodium azide Constituents: BSA, 0.0268% PBS
精製度	Protein G purified
特記事項(精製)	purified from ascites
ポリ/モノ	モノクローナル
クローン名	SPM107
アイソタイプ	IgG1
軽鎖の種類	kappa

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab75697** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
IHC-P		

## 追加情報

IHC-P: 1/500 for 10 min at RT. Perform heat mediated antigen retrieval before commencing with IHC staining protocol, by boiling tissue sections in 10mM citrate, pH 6.0 for 10 min followed by cooling at RT for 20 min.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

## ターゲット情報


### 関連性

Follicle stimulating hormone enables ovarian folliculogenesis to the antral follicle stage and is essential for Sertoli cell proliferation and maintenance of sperm quality in the testis. Members of the pituitary glycoprotein hormone family, of which FSH is one (see also luteinizing hormone, chorionic gonadotropin, and thyroid stimulating hormone), consist of a shared alpha chain and a beta chain encoded by a separate gene. The FSHB gene encodes the beta subunit of follicle stimulating hormone.

### 細胞内局在

Secreted

## 画像



Human anterior pituitary stained with 1/500 ab75697 for 10 minutes at RT.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-FSH antibody [SPM107] (ab75697)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors