

Product datasheet

Anti-ZNF677 antibody ab76054

画像数 1

製品の概要

製品名	Anti-ZNF677 antibody
製品の詳細	Rabbit polyclonal to ZNF677
由来種	Rabbit
アプリケーション	適用あり: WB, ELISA
種交差性	交差種: Human
免疫原	A synthetic peptide corresponding to a region of Human ZNF677
ポジティブ・コントロール	Fetal liver lysate

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: None Constituents: 2% Sucrose, PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab76054** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 µg/ml. Detects a band of approximately 37, 68 kDa (predicted molecular weight: 68 kDa). Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

アプリ ケーション	Abreviews	特記事項
--------------	-----------	------


ELISA
 Use at an assay dependent concentration.
 ELISA titre using peptide based assay 1:312500.

ターゲット情報

関連性 ZNF677 belongs to the krueppel C2H2-type zinc-finger protein family. It contains ten C2H2-type zinc fingers and one KRAB domain. ZNF677 may be involved in transcriptional regulation.

細胞内局在 Nucleus

画像



Western blot - Anti-ZNF677 antibody (ab76054)

Anti-ZNF677 antibody (ab76054) at 1 µg/ml + fetal liver lysate at 10 µg

Secondary

HRP conjugated anti-Rabbit IgG at 1/50000 dilution

Predicted band size: 68 kDa

Observed band size: 37,68 kDa

Gel concentration: 12%

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors