

Product datasheet

Anti-MYO3A antibody ab79982

画像数 1

製品の概要

製品名	Anti-MYO3A antibody
製品の詳細	Rabbit polyclonal to MYO3A
由来種	Rabbit
アプリケーション	適用あり: IHC-P, ELISA
種交差性	交差種: Human
免疫原	KLH conjugated synthetic peptide selected from the C terminal region of human MYO3A.
ポジティブ・コントロール	Human lung carcinoma tissue

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
バッファー	Preservative: 0.09% Sodium azide Constituent: PBS
精製度	Protein A purified
特記事項 (精製)	Purified through a protein A column, eluted with high and low pH buffers and neutralized immediately, followed by dialysis against PBS.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab79982** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
IHC-P		1/50.
ELISA		1/1000.

ターゲット情報

機能	Probable actin-based motor with a protein kinase activity. Probably plays a role in vision and hearing.
組織特異性	Strongest expression in retina, retinal pigment epithelial cells, cochlea and pancreas.
関連疾患	Deafness, autosomal recessive, 30
配列類似性	In the N-terminal section; belongs to the protein kinase superfamily. STE Ser/Thr protein kinase family. Contains 3 IQ domains. Contains 1 myosin head-like domain. Contains 1 protein kinase domain.
細胞内局在	Cytoplasm > cytoskeleton.

画像



ab79982, at a 1/50 dilution, staining MYO3A in formalin fixed, paraffin embedded human lung carcinoma tissue by Immunohistochemistry. ab79982 was peroxidase conjugated to the secondary antibody, followed by DAB staining.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-MYO3A antibody (ab79982)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors