

## Product datasheet

# Goat Anti-Human IgA H&L ab8499

### 製品の概要

製品名	Goat Anti-Human IgA H&L
由来種	Goat
ターゲット生物種	Human
特異性	Antibodies were raised in goat to human IgA. There is low background fluorescence in normal tissue in contrast to cases of autoimmune disease where IgA is deposited. No cross reactivity with albumine fractions of human serum. No cross reactivity with albumine fractions of animal species.
アプリケーション	<b>適用あり:</b> ELISA, Dot blot

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	50% glycerol, 0.1 M NaCl, 0.01 M sodium phosphate at pH 7.5.
精製度	Immunogen affinity purified
特記事項(精製)	From the crude polyclonal the crossreactive antibodies were extracted by incubation with Sepharose-bound human IgG and IgM. Specific antibodies were absorbed by incubation with Sepharose-bound human IgA. Specific antibodies were eluted by acidic buffer at pH 2.5 followed by neutralisation and dialysis. After repeated binding with immobilized human IgA a minimum of 65% protein bound.
ポリ/モノ	ポリクローナル
軽鎖の種類	unknown

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab8499** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
ELISA		

Dot blot

**追加情報**

Recommended dilution as secondary antibody in ELISA is 1:1000.

Recommended dilution as secondary antibody in dot-blotting is 1:100.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

**Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

**Terms and conditions**

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors