

Product datasheet

Anti-KCNK15 antibody ab85511

画像数 1

製品の概要

製品名	Anti-KCNK15 antibody
製品の詳細	Rabbit polyclonal to KCNK15
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Human 交差が予測される動物種: Cat, Dog, Pig 
免疫原	Synthetic peptide corresponding to a region within the internal sequence amino acids 281-330 (ARSVGSASVF CHVHKLERCA RDNLGFSPPS SPGVVRRGGQA PRPGARWKS) of Human KCNK15, NP_071753. Run BLAST with ExPASy Run BLAST with NCBI
ポジティブ・コントロール	Human fetal muscle lysate

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: None Constituents: 2% Sucrose, PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab85511** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 36 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

ターゲット情報

機能	Probable potassium channel subunit. No channel activity observed in heterologous systems. May need to associate with another protein to form a functional channel.
組織特異性	Detected in pancreas, heart, placenta, lung, liver, kidney, ovary, testis, skeletal muscle and adrenal gland, and at lower levels in prostate, spleen and thyroid gland.
配列類似性	Belongs to the two pore domain potassium channel (TC 1.A.1.8) family.
細胞内局在	Membrane.

画像



Anti-KCNK15 antibody (ab85511) at 1 µg/ml
(in 5% skim milk / PBS buffer) + Human fetal
muscle lysate at 10 µg

Secondary

HRP conjugated anti-Rabbit IgG at 1/50000
dilution

Predicted band size: 36 kDa

Observed band size: 39 kDa

Additional bands at: 32 kDa. We are unsure
as to the identity of these extra bands.

Gel concentration 12%

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors