

Product datasheet

Anti-HOXB8 antibody ab85745

画像数 1

製品の概要

製品名	Anti-HOXB8 antibody
製品の詳細	Rabbit polyclonal to HOXB8
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Human 交差が予測される動物種: Dog, Pig
免疫原	Synthetic peptide corresponding to a region within the internal sequence amino acids 194-243 (FQNRRMKWKK ENNKDKFPSS KCEQEELEKQ KLERAPEAAD EGDAQKGDKK) of Human HOXB8, NP_076921. Run BLAST with ExPASy Run BLAST with NCBI
ポジティブ・コントロール	Human fetal muscle lysate

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: None Constituents: 2% Sucrose, PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab85745** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリ ケーション	Abreviews	特記事項
--------------	-----------	------

WB Use a concentration of 1 µg/ml. Detects a band of approximately 28 kDa (predicted molecular weight: 28 kDa). Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

ターゲット情報

機能	Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis.
配列類似性	Belongs to the Antp homeobox family. Contains 1 homeobox DNA-binding domain.
発生段階	Expressed in whole embryos and fetuses at 5-9 weeks from conception.
細胞内局在	Nucleus.

画像



Anti-HOXB8 antibody (ab85745) at 1 µg/ml + fetal muscle lysate at 10 µg

Secondary

HRP conjugated anti-Rabbit IgG at 1/50000 dilution

Predicted band size: 28 kDa

Observed band size: 28 kDa

Additional bands at: 90 kDa. We are unsure as to the identity of these extra bands.

Gel concentration 12%

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors