

## Product datasheet

# Anti-ZNF280D antibody ab85881

### 画像数 1

#### 製品の概要

製品名	Anti-ZNF280D antibody
製品の詳細	Rabbit polyclonal to ZNF280D
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Human 交差が予測される動物種: Rabbit, Horse, Cat, Dog, Pig 
免疫原	Synthetic peptide corresponding to a region within the internal sequence amino acids 539-588 (STLQLSPPRT KNITAKNPAK SNTSKPNTVK SNASKPNTSK PNGSKSKYKP) of Human ZNF280D (NP_060131). <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
ポジティブ・コントロール	Fetal muscle lysate.

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: None Constituents: 2% Sucrose, PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab85881** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリ ケーショ ン	Abreviews	特記事項
------------------	-----------	------

WB Use a concentration of 1 µg/ml. Detects a band of approximately 109 kDa (predicted molecular weight: 109 kDa). Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.


## ターゲット情報

**機能** May function as a transcription factor.

**配列類似性** Contains 5 C2H2-type zinc fingers.

**細胞内局在** Nucleus.

## 画像



Western blot - Anti-ZNF280D antibody (ab85881)

Anti-ZNF280D antibody (ab85881) at 1 µg/ml  
(5% skim milk / PBS buffer) + fetal muscle  
lysate at 10 µg

### Secondary

HRP conjugated anti-Rabbit IgG at 1/50000  
dilution

**Predicted band size:** 109 kDa

**Observed band size:** 109 kDa

Gel concentration 6-18%

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors