

Product datasheet

Anti-TULP2 antibody ab86405

画像数 2

製品の概要

製品名	Anti-TULP2 antibody
製品の詳細	Rabbit polyclonal to TULP2
由来種	Rabbit
アプリケーション	適用あり: IHC-P, WB
種交差性	交差種: Human 交差が予測される動物種: Mouse, Rat, Horse, Cow, Cat
免疫原	Synthetic peptide corresponding to a region within N terminal amino acids 2 - 51 (SQDNDTLMRDILGHELAAMRLQKLEQQRRLEKKQRQKRQELLMVQANPD) of Human TULP2 (NP_003314) Run BLAST with Run BLAST with
ポジティブ・コントロール	HeLa cell lysate

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: 2% Sucrose, PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab86405** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
IHC-P		Use a concentration of 5 µg/ml.
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 59 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

ターゲット情報

組織特異性	Strongly expressed in testis. Also expressed in retina. Expressed in cancer cell lines.
配列類似性	Belongs to the TUB family.
細胞内局在	Cytoplasm. Secreted. Does not have a cleavable signal peptide and is secreted by a non-conventional pathway.

画像




Anti-TULP2 antibody (ab86405) at 1 µg/ml +
HeLa cell lysate at 10 µg

Secondary

HRP conjugated anti-Rabbit IgG at 1/50000
dilution

Predicted band size: 59 kDa

Gel concentration: 12%



Immunohistochemistry (Formalin/PFA-fixed
paraffin-embedded sections) analysis of
human colon tissue labelling TULP2 with
ab86405 at 5.0µg/ml.

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors