


Product datasheet

Anti-SNRPD2 antibody ab86519

★★★★☆ 1 Abreviews [画像数 1](#)

製品の概要

製品名	Anti-SNRPD2 antibody
製品の詳細	Rabbit polyclonal to SNRPD2
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Human 交差が予測される動物種: Mouse, Rat, Rabbit, Goat, Horse, Guinea pig, Cow, Cat, Dog, Drosophila melanogaster, Zebrafish 
免疫原	Synthetic peptide within residues: NVKEMWTEVP KSGKGKKKSK PVNKDRYISK MFLRGDSVIV VLRNPLIAGK, corresponding to C terminal amino acids 69-118 of Human SNRPD2 (NP_004588) Run BLAST with ExPASy Run BLAST with NCBI
ポジティブ・コントロール	HepG2 cell lysate.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: 2% Sucrose, PBS
精製度	Immunogen affinity purified
特記事項(精製)	ab86519 is purified by a peptide affinity chromatography method.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab86519** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリ ケーショ ン	Abreviews	特記事項
------------------	-----------	------

WB ★★★★☆ Use a concentration of 1 µg/ml. Detects a band of approximately 13 kDa (predicted molecular weight: 13 kDa). Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

ターゲット情報

機能 Required for pre-mRNA splicing. Required for snRNP biogenesis.
配列類似性 Belongs to the snRNP core protein family.
細胞内局在 Nucleus.

画像



Anti-SNRPD2 antibody (ab86519) at 1 µg/ml
+ HepG2 cell lysate at 10 µg

Secondary

anti-Rabbit IgG HRP at 1/50000 dilution

Predicted band size: 13 kDa

Observed band size: 13 kDa

Western blot - Anti-SNRPD2 antibody (ab86519)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors