

Product datasheet

Anti-RASL10A antibody ab86981

画像数 1

製品の概要

製品名	Anti-RASL10A antibody
製品の詳細	Rabbit polyclonal to RASL10A
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Human 交差が予測される動物種: Horse 
免疫原	Synthetic peptide corresponding to a region within the N terminal amino acids 36-85 (TDGPRLYRPA VLLDGAVYDL SIRDGDVAGP GSSPGGPPEEW PDAKDWSLQD) of Human RASL10A (NP_006468). Run BLAST with ExPASy Run BLAST with NCBI
ポジティブ・コントロール	Transfected 293T cell lysate.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: None Constituents: 2% Sucrose, PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab86981** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 23 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

ターゲット情報

機能	Potent inhibitor of cellular proliferation.
組織特異性	Expression appears to be strictly limited to the central nervous system.
配列類似性	Belongs to the small GTPase superfamily. Ras family.
翻訳後修飾	Isoprenylation is essential for nucleolar localization, and the proliferation-inhibiting activity of RASL10A.
細胞内局在	Cell membrane. Nucleus > nucleolus. May cycle in and out of the nucleolus in a GTP-dependent manner.

画像



Western blot - Anti-RASL10A antibody (ab86981)

Anti-RASL10A antibody (ab86981) at 1 µg/ml
+ Transfected 293T cell lysate at 10 µg

Secondary

HRP conjugated anti-Rabbit IgG at 1/50000 dilution

Predicted band size: 23 kDa

Observed band size: 23 kDa

Gel concentration: 12%

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors