

Product datasheet

Anti-CSTF3 antibody ab87158

画像数 1

製品の概要

製品名	Anti-CSTF3 antibody
製品の詳細	Rabbit polyclonal to CSTF3
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Human 交差が予測される動物種: Mouse, Rat, Rabbit, Horse, Chicken, Guinea pig, Cow, Cat, Dog, Zebrafish ▲
免疫原	Synthetic peptide corresponding to a region within N-terminal amino acids 108-157 (VRETKGKLPYSYKEKMAQAYDFALDKIGMEIMSYQMVDYINFLKGVEAV G) of human CSTF3 (NP_001317).
ポジティブ・コントロール	721_B Cell lysate

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: 2% Sucrose, PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab87158** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
----------	-----------	------


アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 83 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

ターゲット情報

関連性 CSTF3 (cleavage stimulation factor subunit 3) is required for polyadenylation and 3'-end cleavage of mammalian pre-mRNAs and is part of the CSTF complex. The CSTF complex is composed of CSTF1 (50 kDa subunit), CSTF2 (64 kDa subunit) and CSTF3 (77 kDa subunit). CSTF3 functions as a homodimer and interacts directly with CSTF1 and CSTF2.

細胞内局在 Nuclear

画像



Anti-CSTF3 antibody (ab87158) at 1 µg/ml +
721_B Cell lysate at 10 µg

Secondary

HRP conjugated anti-Rabbit IgG
at 1: 50,000

Predicted band size: 83 kDa

Observed band size: 83 kDa

Additional bands at: 36 kDa. We are unsure as to the identity of these extra bands.

Gel concentration: 12%

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors