

Product datasheet

Anti-Cytohesin 1 antibody ab88092

画像数 3

製品の概要

製品名	Anti-Cytohesin 1 antibody
製品の詳細	Mouse polyclonal to Cytohesin 1
由来種	Mouse
アプリケーション	適用あり: WB
種交差性	交差種: Mouse, Human
免疫原	Full length, Human Cytohesin 1 (NP_004753.1).
ポジティブ・コントロール	Jurkat, NIH/3T3 and transfected 293T cell lysates.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: 1X PBS, pH 7.2
精製度	Protein A purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション


Our [Abpromise guarantee](#) covers the use of **ab88092** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 µg/ml. Detects a band of approximately 46 kDa (predicted molecular weight: 46 kDa).

機能	Promotes guanine-nucleotide exchange on ARF1 and ARF5. Promotes the activation of ARF through replacement of GDP with GTP.
組織特異性	Ubiquitous.
配列類似性	Contains 1 PH domain. Contains 1 SEC7 domain.
細胞内局在	Cell membrane.

画像




Western blot - Anti-Cytohesin 1 antibody (ab88092)

Anti-Cytohesin 1 antibody (ab88092) at 1 µg/ml + Jurkat cell lysate at 50 µg

Predicted band size: 46 kDa

Observed band size: 46 kDa




Western blot - Anti-Cytohesin 1 antibody (ab88092)

Anti-Cytohesin 1 antibody (ab88092) at 1 µg/ml + NIH/3T3 cell lysate at 50 µg

Predicted band size: 46 kDa

Observed band size: 38 kDa

[why is the actual band size different from the predicted?](#)



Western blot - Anti-Cytohesin 1 antibody (ab88092)

All lanes : Anti-Cytohesin 1 antibody (ab88092) at 1 µg/ml

Lane 1 : 293T cell line (Cytohesin 1 transfected lysate)

Lane 2 : 293T cell line (Non-transfected lysate)

Lysates/proteins at 25 µg per lane.

Predicted band size: 46 kDa

Observed band size: 46 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors