

Product datasheet

Anti-MDS028 antibody ab88325

画像数 4

製品の概要

製品名	Anti-MDS028 antibody
製品の詳細	Mouse monoclonal to MDS028
由来種	Mouse
アプリケーション	適用あり: WB, ELISA, ICC/IF
種交差性	交差種: Human
免疫原	Recombinant full length human MDS028 protein (AAH06552) with a proprietary tag.
ポジティブ・コントロール	WB: HL60 cell lysate ICC/IF: HeLa cell

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: 1X PBS, pH 7.2
精製度	Protein A purified
ポリ/モノ	モノクローナル
アイソタイプ	IgG1
軽鎖の種類	kappa

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab88325** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 - 5 µg/ml. Predicted molecular weight: 49 kDa.
ELISA		Use at an assay dependent dilution.


アプリケーション**Abreviews****特記事項**

ICC/IF

Use a concentration of 10 µg/ml.

ターゲット情報**配列類似性**

Contains 2 FG-GAP repeats.


画像

Western blot - Anti-MDS028 antibody (ab88325)

Anti-MDS028 antibody (ab88325) at 5 µg/ml +
tagged, recombinant human MDS028 protein
(the immunogen) at 0.2 µg

Predicted band size: 49 kDa**Observed band size:** 75 kDa


Molecular weight of the tag alone is 26 kDa.



Western blot - Anti-MDS028 antibody (ab88325)


Anti-MDS028 antibody (ab88325) at 5 µg/ml +
HL60 cell lysate at 50 µg

Predicted band size: 49 kDa**Observed band size:** 49 kDa



ab88325, at 10 $\mu\text{g/ml}$, staining MDS028 in HeLa cell by Immunofluorescence.

Immunocytochemistry/ Immunofluorescence - Anti-MDS028 antibody (ab88325)



Detection limit for recombinant proprietary tagged MDS028 is approximately 3 ng/ml as a capture antibody.

ELISA - Anti-MDS028 antibody (ab88325)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors