

## Product datasheet

# Anti-MCART1 antibody ab88767

### 画像数 2

#### 製品の概要

製品名	Anti-MCART1 antibody
製品の詳細	Mouse polyclonal to MCART1
由来種	Mouse
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Full length Human MCART1 (NP_219480.1, 1 a.a. - 297 a.a)
ポジティブ・コントロール	Human kidney tissue lysate. MCART1 transfected 293T cell lysate.

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: 1X PBS, pH 7.2
精製度	Protein A purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab88767** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 34 kDa.

#### ターゲット情報


配列類似性 Belongs to the mitochondrial carrier family.

Contains 3 Solcar repeats.

細胞内局在

Mitochondrion inner membrane.

画像




Western blot - Anti-MCART1 antibody (ab88767)

Anti-MCART1 antibody (ab88767) at 1 µg/ml + Human kidney lysate at 50 µg

**Predicted band size:** 34 kDa

**Observed band size:** 36 kDa

[why is the actual band size different from the predicted?](#)



Western blot - Anti-MCART1 antibody (ab88767)

**All lanes :** Anti-MCART1 antibody (ab88767) at 1 µg/ml

**Lane 1 :** MCART1 transfected 293T cell lysate

**Lane 2 :** Non-transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

**Predicted band size:** 34 kDa

**Observed band size:** 32 kDa [why is the actual band size different from the predicted?](#)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors

