

Product datasheet

Anti-C9orf156 antibody ab88810

画像数 2

製品の概要

製品名	Anti-C9orf156 antibody
製品の詳細	Mouse polyclonal to C9orf156
由来種	Mouse
アプリケーション	適用あり: WB
種交差性	交差種: Rat, Human
免疫原	Full length Human C9orf156 (AAH02863)
ポジティブ・コントロール	Rat brain lysate. C9orf156 transfected 293T cell lysate.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: PBS, pH 7.2
精製度	Protein A purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab88810** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
WB		1/500 - 1/1000. Predicted molecular weight: 49 kDa.

ターゲット情報

関連性

C9orf156 hydrolyzes acyl-CoA thioesters (in vitro). It has a preference for substrates with medium chain length (C10-C14) and is inactive towards substrates with C18 or C20 aliphatic chains. Its physiological function is not known. C9orf156 interacts with Nef from human immunodeficiency virus (HIV-1).

画像




Western blot - Anti-C9orf156 antibody (ab88810)

Anti-C9orf156 antibody (ab88810) at 1/500 dilution + Rat brain lysate at 50 μ g

Predicted band size: 49 kDa

Observed band size: 49 kDa



Western blot - Anti-C9orf156 antibody (ab88810)

All lanes : Anti-C9orf156 antibody (ab88810) at 1/500 dilution

Lane 1 : C9orf156 transfected 293T cell lysate

Lane 2 : Non-transfected 293T cell lysate

Lysates/proteins at 25 μ g per lane.

Predicted band size: 49 kDa

Observed band size: 52 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors