

## Product datasheet

# Anti-Cyclin T2 antibody ab89450

### 画像数 2

#### 製品の概要

製品名	Anti-Cyclin T2 antibody
製品の詳細	Mouse monoclonal to Cyclin T2
由来種	Mouse
アプリケーション	<b>適用あり:</b> WB, ELISA
種交差性	<b>交差種:</b> Human
免疫原	Recombinant fragment, corresponding to amino acids 264-371 of Human Cyclin T2 (NP_490595)
ポジティブ・コントロール	Jurkat cell lysate.

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: 1X PBS, pH 7.2
精製度	Protein A purified
ポリ/モノ	モノクローナル
アイソタイプ	IgG2a
軽鎖の種類	kappa

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab89450** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 81 kDa.
ELISA		Use at an assay dependent dilution.

## ターゲット情報

機能	Regulatory subunit of the cyclin-dependent kinase pair (CDK9/cyclin T) complex, also called positive transcription elongation factor B (P-TEFB), which is proposed to facilitate the transition from abortive to production elongation by phosphorylating the CTD (carboxy-terminal domain) of the large subunit of RNA polymerase II (RNAP II).
組織特異性	Ubiquitously expressed.
配列類似性	Belongs to the cyclin family. Cyclin C subfamily.
細胞内局在	Nucleus.

## 画像




Western blot - Anti-Cyclin T2 antibody (ab89450)

Anti-Cyclin T2 antibody (ab89450) at 1 µg/ml  
+ Cyclin T2 recombinant fragment at 0.2 µg

**Predicted band size:** 81 kDa

**Observed band size:** 38 kDa

The observed band size may not correspond to the predicted protein molecular weight as the immunogen (recombinant fragment) was used for the test lane.



Western blot - Anti-Cyclin T2 antibody (ab89450)

Anti-Cyclin T2 antibody (ab89450) at 5 µg/ml  
+ Jurkat cell lysate at 50 µg

**Predicted band size:** 81 kDa

**Observed band size:** 74 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors