

Product datasheet

Anti-TSHZ1 antibody ab91230

画像数 1

製品の概要

製品名	Anti-TSHZ1 antibody
製品の詳細	Rabbit polyclonal to TSHZ1
由来種	Rabbit
アプリケーション	適用あり: WB, ELISA
種交差性	交差種: Mouse, Human
免疫原	A KLH conjugated synthetic peptide selected from the Center region of Human TSHZ1.
ポジティブ・コントロール	Mouse kidney tissue lysates

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: 0.09% Sodium Azide Constituents: PBS
精製度	Ammonium Sulphate Precipitation
特記事項 (精製)	ab91230 is purified by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab91230** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		
ELISA		

追加情報

ELISA: 1/1000.

WB: 1/50 - 1/100. Predicted molecular weight: 118 kDa.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

ターゲット情報


関連性

This gene (TSHZ1) encodes a colon cancer antigen which was defined by serological analysis of recombinant cDNA expression libraries. The encoded protein is a member of the teashirt C2H2-type zinc-finger protein family and may be involved in transcriptional regulation of developmental processes.

細胞内局在

Nuclear

画像



Anti-TSHZ1 antibody (ab91230) at 1/50

dilution + mouse kidney tissue lysates at 35 µg

Predicted band size: 118 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors