

## Product datasheet

# Anti-CD300e antibody [UP-H2] (Allophycocyanin) ab91347

### 製品の概要

|          |  |
|----------|--|
| 製品名      | Anti-CD300e antibody [UP-H2] (Allophycocyanin)                 |
| 製品の詳細    | Mouse monoclonal [UP-H2] to CD300e (Allophycocyanin)           |
| 由来種      | Mouse  |
| 標識       | Allophycocyanin. Ex: 645nm, Em: 660nm                          |
| アプリケーション | <b>適用あり:</b> Flow Cyt  |
| 種交差性     | <b>交差種:</b> Human  |
| 免疫原      | Tissue/ cell preparation (Human): IREM-1-HA transfected cells. |

### 製品の特性

|        |  |
|--------|--|
| 製品の状態  | Liquid   |
| 保存方法   | Shipped at 4°C. Store at +4°C.                           |
| バッファー  | Preservative: 0.09% Sodium Azide<br>Constituents: 1% BSA |
| 精製度    | Immunogen affinity purified                              |
| ポリ/モノ  | モノクローナル  |
| クローン名  | UP-H2  |
| アイソタイプ | IgG1   |

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab91347** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| アプリケーション | Abreviews | 特記事項                                |
|----------|-----------|-------------------------------------|
| Flow Cyt |           | Use at an assay dependent dilution. |

### ターゲット情報

|       |   |
|-------|---|
| 機能    | Probably acts as an activating receptor.  |
| 組織特異性 | Present on the surface of mature hematopoietic cells of the monocyte and myeloid lineages (at protein level). |
| 配列類似性 | Belongs to the CD300 family.<br>Contains 1 Ig-like V-type (immunoglobulin-like) domain.                       |
| 翻訳後修飾 | N-glycosylated.   |
| 細胞内局在 | Cell membrane.  |

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors