

## Product datasheet

# Anti-CSGALNACT2 antibody ab93597

### 画像数 1

#### 製品の概要

製品名	Anti-CSGALNACT2 antibody
製品の詳細	Rabbit polyclonal to CSGALNACT2
由来種	Rabbit
特異性	Recognizes isoform 2
アプリケーション	<b>適用あり:</b> WB, ELISA
種交差性	<b>交差種:</b> Human
免疫原	Synthetic peptide selected from the center region of Human CSGALNACT2 conjugated to KLH (NP_061060.3).
ポジティブ・コントロール	293 cell line lysates.

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: 0.09% Sodium Azide Constituents: PBS
精製度	Immunogen affinity purified
特記事項(精製)	ab93597 is purified through a protein A column, followed by peptide affinity purification.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab93597** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		

ELISA

## 追加情報

ELISA: 1/1000.

WB: 1/100 - 1/500. Predicted molecular weight: 38 kDa.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

## ターゲット情報

## 機能

Transfers 1,4-N-acetylgalactosamine (GalNAc) from UDP-GalNAc to the non-reducing end of glucuronic acid (GlcUA). Required for addition of the first GalNAc to the core tetrasaccharide linker and for elongation of chondroitin chains.

## 組織特異性

Ubiquitous.


## 配列類似性

Belongs to the chondroitin N-acetylgalactosaminyltransferase family.

## 細胞内局在

Golgi apparatus &gt; Golgi stack membrane.

## 画像



Anti-CSGALNACT2 antibody (ab93597) at  
1/100 dilution + 293 cell line lysates at 35 µg

**Predicted band size: 38 kDa**

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors