

## Product datasheet

# Anti-CD28 antibody [JJ319] (Biotin) ab95561

### 画像数 1

#### 製品の概要

---

製品名	Anti-CD28 antibody [JJ319] (Biotin)
製品の詳細	Mouse monoclonal [JJ319] to CD28 (Biotin)
由来種	Mouse
標識	Biotin
アプリケーション	<b>適用あり:</b> Flow Cyt
種交差性	<b>交差種:</b> Rat
免疫原	Rat CD28 transfected Mouse A20JB lymphoma
ポジティブ・コントロール	Rat splenocytes

#### 製品の特性

---

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.09% Sodium Azide Constituents: PBS, pH 7.2
精製度	Protein G purified
ポリモノ	モノクローナル
クローン名	JJ319
アイソタイプ	IgG1

#### アプリケーション

---

Our [Abpromise guarantee](#) covers the use of **ab95561** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
----------	-----------	------


Flow Cyt Use 1µg for 10<sup>5</sup> cells.  
 Use at less than or equal to 1 µg per test. A test is defined as the amount (µg) of antibody that will stain a cell sample in a final volume of 100 µL. Cell number should be determined empirically but can range from 105 to 108 cells/test.

[ab18434](#) - Mouse monoclonal IgG1, is suitable for use as an isotype control with this antibody.

### ターゲット情報

機能	Involved in T-cell activation, the induction of cell proliferation and cytokine production and promotion of T-cell survival.
組織特異性	Expressed in T-cells and plasma cells, but not in less mature B-cells.
配列類似性	Contains 1 Ig-like V-type (immunoglobulin-like) domain.
細胞内局在	Membrane.

### 画像



Flow cytometry staining of rat splenocytes with Anti-Rat CD5 FITC and 0.25 µg of Mouse IgG1 ? Isotype Control Biotin (left) or 0.5 µg of ab95561 (right) followed by Streptavidin PE. Total viable cells were used for analysis.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors