

## Product datasheet

# Anti-HERC6 antibody ab99173

### 画像数 1

#### 製品の概要

|          |   |
|----------|---|
| 製品名      | Anti-HERC6 antibody   |
| 製品の詳細    | Rabbit polyclonal to HERC6  |
| 由来種      | Rabbit  |
| アプリケーション | 適用あり: WB  |
| 種交差性     | 交差種: Human<br>交差が予測される動物種: Mouse, Rat, Rabbit, Horse, Chicken, Guinea pig, Cow, Cat, Dog, Pig |



|     |  |
|-----|--|
| 免疫原 | Synthetic peptide corresponding to a region within N terminal amino acids 143-192 (LSKDSQVFSW GKNSHGQLGL GKEFPSQASP QRVRSLEGIP LAQVAAGGAH) of Human HERC6 (NP_060382). |
|-----|--|

Run BLAST with **ExPASy**
 Run BLAST with **NCBI**

#### 製品の特性

|        |  |
|--------|--|
| 製品の状態  | Liquid   |
| 保存方法   | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles. |
| バッファー  | Preservative: None<br>Constituents: 2% Sucrose, PBS  |
| 精製度    | Immunogen affinity purified  |
| ポリ/モノ  | ポリクローナル  |
| アイソタイプ | IgG  |

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab99173** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|------|
|          |           |      |

| アプリケーション | Abreviews | 特記事項  |
|----------|-----------|---|
| WB       |           | Use a concentration of 1 µg/ml. Predicted molecular weight: 115 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T. |

## ターゲット情報

|       |  |
|-------|--|
| 機能    | E3 ubiquitin-protein ligase which accepts ubiquitin from an E2 ubiquitin-conjugating enzyme in the form of a thioester and then directly transfers the ubiquitin to targeted substrates. |
| 組織特異性 | Detected in brain, heart, placenta and testis.   |
| パスウェイ | Protein modification; protein ubiquitination.  |
| 配列類似性 | Contains 1 HECT (E6AP-type E3 ubiquitin-protein ligase) domain.<br>Contains 5 RCC1 repeats.  |
| 発生段階  | Expressed in fetal brain.  |
| 細胞内局在 | Cytoplasm > cytosol.   |

## 画像



Western blot analysis of ACHN cell lysate labeling HERC6 with ab99173 at 1µg/ml.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors