

## Product datasheet

# Mouse monoclonal [G1 1C5] Anti-Rat IgG1 heavy chain (FITC) ab99659

### 製品の概要

|          |   |
|----------|---|
| 製品名      | Mouse monoclonal [G1 1C5] Anti-Rat IgG1 heavy chain (FITC)        |
| 由来種      | Mouse   |
| ターゲット生物種 | Rat   |
| 特異性      | Reacts with the gamma 1 heavy chain (Fc) of Rat IgG1.             |
| アプリケーション | <b>適用あり:</b> ICC/IF, Flow Cyt, ELISA                              |
| 免疫原      | The details of the immunogen for this antibody are not available. |
| 標識       | FITC. Ex: 493nm, Em: 528nm  |

### 製品の特性

|        |  |
|--------|--|
| 製品の状態  | Liquid   |
| 保存方法   | Shipped at 4°C. Store at +4°C.                       |
| バッファー  | Preservative: 0.1% Sodium Azide<br>Constituents: PBS |
| 精製度    | Protein G purified                                   |
| ポリ/モノ  | モノクローナル  |
| クローン名  | G1 1C5   |
| アイソタイプ | IgG1   |
| 軽鎖の種類  | kappa  |

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab99659** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| アプリケーション | Abreviews | 特記事項                                     |
|----------|-----------|--|
| ICC/IF   |           | Use at an assay dependent concentration. |

| アプリケーション | Abreviews | 特記事項  |
|----------|-----------|---|
| Flow Cyt |           | Use 0.3µg for 10 <sup>6</sup> cells.<br><a href="#">ab106163</a> - Mouse monoclonal IgG1, is suitable for use as an isotype control with this antibody. |
| ELISA    |           | Use at an assay dependent concentration.  |

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors