

Product datasheet

Mouse monoclonal [2A9] Anti-Rabbit IgG heavy chain (Alkaline Phosphatase) ab99701

製品の概要

製品名	Mouse monoclonal [2A9] Anti-Rabbit IgG heavy chain (Alkaline Phosphatase)
由来種	Mouse
ターゲット生物種	Rabbit
特異性	ab99701 reacts with the heavy chain of all subclasses of Rabbit IgG.
アプリケーション	適用あり: ELISA, Flow Cyt
標識	Alkaline Phosphatase

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.1% Sodium Azide Constituents: 50% Glycerol, 50mM Tris, 1mM Magnesium chloride, pH 8.0
精製度	Protein G purified
ポリ/モノ	モノクローナル
クローン名	2A9
アイソタイプ	IgG1

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab99701** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
ELISA		1/2000 - 1/4000.
Flow Cyt		Use at an assay dependent dilution.

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors