

## Product datasheet

# alpha 1 Adrenergic Receptor peptide ab41797

## 1 References

### 製品の概要

製品名 alpha 1 Adrenergic Receptor peptide

### 製品の詳細

由来 Synthetic

### アミノ酸配列

配列 K-FSREKKAAT

### 特性

Our [Abpromise guarantee](#) covers the use of **ab41797** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション Neutralising  
Blocking - Blocking peptide for Anti-alpha 1 Adrenergic Receptor antibody ([ab3462](#))

製品の状態 Liquid

### 前処理および保存

保存方法および安定性 Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

### 関連情報

機能 This alpha-adrenergic receptor mediates its effect through the influx of extracellular calcium.

配列類似性 Belongs to the G-protein coupled receptor 1 family. Adrenergic receptor subfamily. ADRA1D sub-subfamily.

細胞内局在 Cell membrane.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors