

Product datasheet

Anti-beta Amyloid 37-42 antibody ab34376

2 References [画像数 1](#)

製品の概要

製品名	Anti-beta Amyloid 37-42 antibody
製品の詳細	Rabbit polyclonal to beta Amyloid 37-42
由来種	Rabbit
アプリケーション	適用あり: ELISA, WB
種交差性	交差種: Human
免疫原	Synthetic peptide corresponding to beta Amyloid 37-42 (C terminal).

製品の特性


製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
バッファー	Preservative: None Constituents: PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab34376** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
ELISA		1/3000.
WB		1/400. Detects a band of approximately 85 kDa (predicted molecular weight: 87 kDa). 1/1000 for analysis of the recombinant protein, 1/400 for tissue extracts or cell lysates.



Western blot - Anti-beta Amyloid 37-42 antibody (ab34376)

Anti-beta Amyloid 37-42 antibody (ab34376)
at 1 µg/ml + Human Plasma Total Protein
Lysate at 10 µg

Secondary

Goat Anti-Rabbit IgG H&L (HRP) preadsorbed (ab97080) at 1/5000 dilution

Developed using the ECL technique.

Performed under reducing conditions.

Predicted band size: 87 kDa

Observed band size: 85 kDa

Additional bands at: 28 kDa, 55 kDa. We are unsure as to the identity of these extra bands.

Exposure time: 90 seconds

The band observed at 85 kDa could potentially be a cleaved form of beta Amyloid 37-42 due to the presence of a 17 amino acid signal peptide.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors