

## Product datasheet

# Donkey Anti-Sheep IgG H&L ab6895

### 1 References

#### 製品の概要

製品名	Donkey Anti-Sheep IgG H&L
由来種	Donkey
ターゲット生物種	Sheep
アプリケーション	<b>適用あり:</b> WB, ELISA, ICC/IF, IHC-P, IHC-Fr, IP, Immunodiffusion, Conjugation
免疫原	Full length native Sheep IgG (purified).

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.01% Sodium azide Constituents: 0.42% Potassium phosphate, 0.87% Sodium chloride
精製度	Affinity purified
特記事項 (標識)	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Sheep IgG coupled to agarose
ポリ/モノ	ポリクローナル

#### アプリケーション

**The Abpromise guarantee** [Abpromise保証は](#)、次のテスト済みアプリケーションにおけるab6895の使用に適用されます

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		
ELISA		
ICC/IF		

アプリケーション	Abreviews	特記事項
IHC-P		
IHC-Fr		
IP		
Immunodiffusion		
Conjugation		

#### 追加情報

Conjugation: Use at an assay dependant dilution.

ID: Use at an assay dependant dilution.

IHC-P: Use at an assay dependant dilution.

IHC-Fr: Use at an assay dependant dilution.

IP: Use at an assay dependant dilution.

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

#### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors