

## Product datasheet

# Goat Anti-Armenian hamster IgG H&L (FITC) ab5739

## 4 References

### 製品の概要

---

製品名	Goat Anti-Armenian hamster IgG H&L (FITC)
由来種	Goat
ターゲット生物種	Armenian hamster
特異性	Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Fluorescein, anti-Goat Serum, Armenian Hamster IgG and Armenian Hamster Serum.
アプリケーション	<b>適用あり:</b> ICC/IF, IHC-P, ELISA, Flow Cyt, Immunomicroscopy
免疫原	Armenian Hamster IgG whole molecule
標識	FITC. Ex: 493nm, Em: 528nm

### 製品の特性

---

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Preservative: 0.01% Sodium azide Constituents: 0.42% Potassium phosphate, 0.87% Sodium chloride, 1% BSA  BSA is protease free
特記事項 (精製)	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Armenian Hamster IgG coupled to agarose beads.
特記事項 (標識)	Fluorescein isothiocyanate (FITC) (Molecular Weight 390 daltons). Absorption Wavelength: 495 nm Emission Wavelength: 528 nm Fluorochrome/Protein Ratio: 4.3 moles FITC per mole of Goat IgG
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション

---

**The Abpromise guarantee** [Abpromise保証は](#)、次のテスト済みアプリケーションにおけるab5739の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ICC/IF		1/1000 - 1/5000.
IHC-P		Use at an assay dependent concentration. PubMed: 24146823
ELISA		1/10000 - 1/50000.
Flow Cyt		1/500 - 1/2500.
Immunomicroscopy		Use at an assay dependent concentration.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors