

# Mesenchymal Stromal Cell Marker (CD44, CD45, CD90 CD29, CD105) Antibody Panel - Human ab93758

[10 References](#) [画像数 2](#)

### 製品の概要

---

**製品名** Mesenchymal Stromal Cell Marker (CD44, CD45, CD90 CD29, CD105) Antibody Panel - Human

**製品の概要** Mesenchymal Stromal Cell Marker (CD44, CD45, CD90 CD29, CD105) Antibody Panel ab93758 contains multiple trial-sized versions of clones against key human mesenchymal stromal cell markers, and includes a CD44 mouse monoclonal, CD45 rabbit polyclonal, CD90 mouse monoclonal, CD29 rabbit monoclonal and CD105 mouse monoclonal antibody.

For guidelines on how to use each antibody within the panel, please consult the individual datasheet for each antibody.

Panel contains:

Anti-CD44 antibody [F10-44-2] ([ab61124](#))

Anti-CD45 antibody - Hematopoietic Stem Cell Marker ([ab10559](#))

Anti-CD90 / Thy1 antibody [AF-9] ([ab23894](#))

Anti-CD105 antibody [105C02] ([ab44967](#))

Recombinant Anti-Integrin beta 1 antibody [EP1041Y] ([ab52971](#))

The Human Multipotent Mesenchymal Stromal Cell Marker Antibody Panel is designed for the characterization of cultured or isolated human multipotent mesenchymal stromal cells. The panel contains a group of antibodies for the positive (CD105, CD29, CD44, CD90) and negative (CD45) selection of this lineage. Multipotent mesenchymal stromal cells (MSCs) represent a heterogeneous subset of stromal cells. They can be bone marrow-derived or isolated from many other different adult tissues including periosteum, trabecular bone, adipose tissue, and synovium. They exhibit the potential to give rise to cells of diverse lineages including adipocytes, chondrocytes and osteocytes.

### 特記事項

**Explore our range of antibody sample panels** designed to provide you with a variety of trial-size antibodies in a convenient and cost-effective format.


## 製品の特性

### 保存方法

Store at -20°C. Please refer to protocols.

内容	1 units
<b><u>ab44967 - Anti-CD105 antibody [105C02]</u></b>	1 x 50µl
<b><u>ab6124 - Anti-CD44 antibody [F10-44-2]</u></b>	1 x 50µg
<b><u>ab10559 - Anti-CD45 antibody</u></b>	1 x 50µg
<b><u>ab23894 - Anti-CD90 / Thy1 antibody [AF-9]</u></b>	1 x 50µg
<b><u>ab134179 - Anti-Integrin beta 1 antibody [EPR1040Y]</u></b>	1 x 10µl

## 画像



Ab6124 staining CD44 on Human Tissue sections (Skin tumor) by IHC-P. Sections were Formaldehyde fixed and subjected to heat mediated antigen retrieval in citric acid prior to blocking with 10% serum for 30 mins at 25°C. The primary antibody was diluted 1/100 in PBS and incubated with the sample overnight at 4°C. A TRITC non-Abcam secondary antibody was used.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Human Mesenchymal Stromal Cell Marker Panel (CD44, CD45, CD90 CD29 and CD105) (ab93758)

### Why choose a recombinant antibody?



**Research with confidence**  
Consistent and reproducible results



**Long-term and scalable supply**  
Recombinant technology



**Success from the first experiment**  
Confirmed specificity



**Ethical standards compliant**  
Animal-free production

Mesenchymal Stromal Cell Marker (CD44, CD45, CD90 CD29, CD105) Antibody Panel - Human (ab93758)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors