

96-Well Microplate Pack ab210903

1 References

製品の概要

製品名	96-Well Microplate Pack
製品の概要	<p>96-Well Microplate Pack is recommended for use with our Matched Antibody Pair kits. The 96-Well plate pack contains 10 solid microplates (not strip well) which are uncoated, high binding microplates with 10 plate seals. The plate can be stored at 4°C or at room temperature.</p> <p>The 96-Well Plate Pack is also available as part of our Matched Antibody Pair Kit Accessory Pack (ab210905).</p>

製品の特性

保存方法 Store at +4°C. Please refer to protocols.

内容	10 packs
96-Well Microplate	10 x 96 tests
plate seal	10 units

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors