

Eph receptor A3 overexpression 293T lysate (whole cell) ab94111

画像数 2

医薬用外劇物

製品の概要

製品名	Eph receptor A3 overexpression 293T lysate (whole cell)
特記事項	ab94111 is a 293T cell transfected lysate in which Human Eph receptor A3 has been transiently over-expressed using a pCMV-Eph receptor A3 plasmid. The lysate is provided in 1X Sample Buffer.
アプリケーション	適用あり: WB

製品の特性


Mycoplasma free	Yes
製品の状態	Liquid
保存方法	Shipped on dry ice. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Constituents: 0.01% Bromophenol blue, 2.3% Beta mercaptoethanol, 2% Sodium lauryl sulfate, 0.788% Tris HCl, 10% Glycerol (glycerin, glycerine)
背景	Function: Receptor for members of the ephrin-A family. Binds to ephrin-A2, -A3, -A4 and -A5. Could play a role in lymphoid function. Tissue specificity: Widely expressed. Highest level in placenta. Disease: Defects in EPHA3 may be a cause of colorectal cancer (CRC) [MIM:114500]. Similarity: Belongs to the protein kinase superfamily. Tyr protein kinase family. Ephrin receptor subfamily. Contains 2 fibronectin type-III domains. Contains 1 protein kinase domain. Contains 1 SAM (sterile alpha motif) domain.

アプリケーション

The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab94111の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。


アプリケーション	Abreviews	特記事項
WB		Use at an assay dependent dilution.

画像



SDS-PAGE - Eph receptor A3 overexpression 293T lysate (whole cell) (ab94111)

ab94111 at 15µg/lane on an SDS-PAGE gel



Western blot - Eph receptor A3 overexpression 293T lysate (whole cell) (ab94111)

All lanes : Anti-Eph receptor A3 antibody (**ab54623**) at 1/500 dilution

Lane 1 : Eph receptor A3 overexpression 293T lysate (whole cell) (ab94111)

Lane 2 : 293T non Transfected Lysate

Lysates/proteins at 25 µg per lane.

Secondary

All lanes : Goat Anti-mouse IgG (H and L) HRP conjugated at 1/2500 dilution

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors

